


## Accuray Precision® Plan Overview

|  |  |  |                               |
|--|--|--|-------------------------------|
| <b>Patient Name:</b>                   | FASTRE, RENE                                 | <b>Date Of Birth:</b>                      | 14 Apr 1931                   |
| <b>Gender:</b>                         | Male   | <b>Patient Key:</b>                        | 409536                        |
| <b>Medical ID:</b>                     | 202305822                                    | <b>Version:</b>                            | Accuray Precision 3.3.1.3 [2] |
| <b>Plan Name:</b>                      | CanalAnal                                    | <b>Prescription:</b>                       | Median of PTV_HDOPT, 55.00 Gy |
| <b>Plan Status:</b>                    | Approved                                     | <b>Prescribed Dose per Fraction (Gy):</b>  | 1.83                          |
| <b>Plan Type:</b>                      | Standard                                     | <b>Planned Fractions:</b>                  | 30                            |
| <b>Treatment Machine / Serial No.:</b> | 4010012_1000MU 1180 MU/min, rev 30 / 4010012 | <b>Maximum Dose (Gy):</b>                  | 57.25                         |
| <b>Plan Saved Date:</b>                | 23 Jun 2023, 3:11:22 PM (hr:min:sec)         | <b>MU Planned / MU per Fraction:</b>       | 75531.8 / 2517.7              |
| <b>Field Width (cm):</b>               | 5.0, Dynamic                                 | <b>Red Lasers Offset (IECf, mm):</b>       | X: 7.76 Y: 18.63 Z: -15.44    |
| <b>Pitch:</b>                          | 0.430  | <b>Couch Travel (mm):</b>                  | 267.7                         |
| <b>Modulation Factor / Actual:</b>     | 2.500 / 1.5                                  | <b>Couch Speed (mm/sec):</b>               | 1.8                           |
| <b>Sinogram Segments:</b>              | 4.1  | <b>Gantry Period (sec):</b>                | 11.8                          |
| <b>Beam On Time (sec):</b>             | 145.5  | <b>Gantry Rotations:</b>                   | 12.2                          |
| <b>Plan Delivery Type / Mode:</b>      | Helical / IMRT                               | <b>Dose Calculation Resolution / Type:</b> | High / Final                  |
| <b>Dose Calculation Algorithm:</b>     | Convolution-Superposition                    | <b>Spacing (IECp, mm)</b>                  | X: 0.98 Y: 2.50 Z: 0.98       |
| <b>Planning Method:</b>                | Classic                                      | <b>Size (IECp, voxel)</b>                  | X: 559 Y: 203 Z: 512          |
| <b>Tracking Method:</b>                | None   | <b>Reference Point (IECp, mm):</b>         | X: 45.9 Y: -130.0 Z: 6.8      |
| <b>Image Angles (deg):</b>             |  | <b>Reference Dose (Gy):</b>                | 57.25                         |
| <b>Tracking Target VOI:</b>            |  |  |                               |

DVH: Absolute Dose



## VOI

### VOI List

| VOI                  | Volume (cm³) | Min (Gy) | Mean (Gy) | Max (Gy) | CI   | nCI  | HI   | Coverage % | Beam Inter. |
|----------------------|--------------|----------|-----------|----------|------|------|------|------------|-------------|
| Grele                | 1587.01      | 0.50     | 5.38      | 42.00    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| QueueDeCheval        | 36.14        | 0.38     | 3.80      | 34.90    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Rectum               | 289.64       | 37.33    | 51.18     | 56.32    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Sigmoide             | 32.75        | 1.20     | 15.95     | 41.85    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| TeteFemorDte         | 230.91       | 12.12    | 20.87     | 35.78    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| TeteFemorGche        | 223.09       | 11.58    | 20.89     | 34.89    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Vessie               | 1218.77      | 1.04     | 22.27     | 46.65    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| BODY                 | 35823.80     | 0.10     | 10.69     | 57.25    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Sacro iliaques       | 387.70       | 1.54     | 17.86     | 42.50    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| OGE                  | 372.40       | 2.20     | 21.97     | 52.63    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| CTV BD               | 1330.32      | 31.52    | 44.93     | 57.25    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| PTV BD               | 2096.97      | 25.86    | 43.58     | 57.25    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| PTV HD               | 424.61       | 48.22    | 54.81     | 57.25    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| GTVp                 | 126.77       | 53.90    | 55.02     | 56.30    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| CTV HD OK            | 290.25       | 53.90    | 55.22     | 57.25    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Grele-PTVs           | 1557.92      | 0.50     | 4.74      | 41.17    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| OGE-PTVs             | 371.33       | 2.20     | 21.93     | 51.86    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| PTV_BDOPT            | 2094.36      | 25.86    | 43.58     | 57.25    | 1.00 | 9.22 | 1.04 | 10.85      | n/a         |
| PTV_HDOPT            | 423.93       | 48.22    | 54.81     | 57.25    | 1.00 | 1.87 | 1.04 | 53.59      | n/a         |
| QueueDeCheval+5      | 113.43       | 0.35     | 5.78      | 39.94    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| RingPTVBD            | 1834.64      | 7.90     | 34.98     | 54.82    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| RingPTVHD            | 372.87       | 20.82    | 44.27     | 54.82    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Sigmoide-PTVs        | 27.18        | 1.20     | 11.39     | 40.87    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| Vessie-PTVs          | 1144.45      | 1.04     | 21.16     | 46.65    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| CHAUFFE              | 852.11       | 13.49    | 26.22     | 39.53    | n/a  | n/a  | n/a  | n/a        | Allowed     |
| All Target Regions   | n/a          | 25.86    | 43.58     | 57.25    | 0.00 | 0.00 | 0.00 | 0.00       | n/a         |
| All Critical Regions | n/a          | 0.10     | 10.71     | 57.25    | 0.00 | 0.00 | 0.00 | 0.00       | n/a         |
| Soft Tissue          | n/a          | 0.10     | 2.49      | 32.49    | 0.00 | 0.00 | 0.00 | 0.00       | n/a         |

### Dx Vx Values

| DVH           | Dose (Gy) | Dose (% of 55.00 Gy Rx Dose) | Volume (cm³) | Volume (%) | Criteria                      | Value  |   |
|---------------|-----------|------------------------------|--------------|------------|-------------------------------|--------|---|
| PTV_HDOPT     | 52.40     | 95.3                         | 402.73       | 95.0       | Dose (Gy) >= 57.00            | 52.40  | ✗ |
| PTV_BDOPT     | 38.81     | 70.6                         | 1989.64      | 95.0       | Dose (Gy) >= 45.60            | 38.81  | ✗ |
| PTV_HDOPT     | 56.15     | 102.1                        | 8.48         | 2.0        | Dose (Gy) <= 64.20            | 56.15  | ✓ |
| Sigmoide      | 40.38     | 73.4                         | 1.00         | 3.1        | Dose (Gy) <= 60.00            | 40.38  | ✓ |
| Sigmoide      | 40.00     | 72.7                         | 2.01         | 6.1        | Volume (%) <= 40.0, Max: 60.0 | 6.1    | ✓ |
| Vessie        | 41.29     | 75.1                         | 1.00         | 0.1        | Dose (Gy) <= 60.00            | 41.29  | ✓ |
| Vessie        | 40.00     | 72.7                         | 17.90        | 1.5        | Volume (%) <= 40.0, Max: 60.0 | 1.5    | ✓ |
| Vessie        | 60.00     | 109.1                        | 0.00         | 0.0        | Volume (%) <= 50.0            | 0.0    | ✓ |
| TeteFemorDte  | 50.00     | 90.9                         | 0.00         | 0.0        | Volume (%) <= 10.0            | 0.0    | ✓ |
| TeteFemorGche | 50.00     | 90.9                         | 0.00         | 0.0        | Volume (%) <= 10.0            | 0.0    | ✓ |
| Grele         | 42.00     | 76.4                         | 0.00         | 0.0        | Dose (Gy) <= 59.00            | 42.00  | ✓ |
| Grele         | 15.00     | 27.3                         | 212.40       | 13.4       | Volume (cm³) <= 830.00        | 212.40 | ✓ |
| Grele         | 45.00     | 81.8                         | 0.00         | 0.0        | Volume (cm³) <= 150.00        | 0.00   | ✓ |
| QueueDeCheval | 34.90     | 63.5                         | 0.00         | 0.0        | Dose (Gy) <= 54.00            | 34.90  | ✓ |

Density

Density Model: IUC-120kV

| HU    | Relative Electron Density | Mass Density (g/cm³) |
|-------|---------------------------|----------------------|
| -939  | *                         | 0.001                |
| -672  | *                         | 0.290                |
| -509  | *                         | 0.480                |
| 13    | *                         | 1.000                |
| 226   | *                         | 1.146                |
| 455   | *                         | 1.331                |
| 837   | *                         | 1.558                |
| 1214  | *                         | 1.822                |
| 30884 | *                         | 22.599               |

## Approval and ICD

### Approved by:

**Name:** MODESTO, Anouchka  
**User ID:** modesto  
**Approval Date:** 23 Jun 2023, 3:11:51 PM (hr:min:sec)

### Specialists

**Radiation Oncologist:**  
**Surgeon:**  
**Referring Physician:**  
**Physicist:**





**Treatment Site:** IUC-TOULOUSE  
**Treatment Site ICD:** IUC-TOULOUSE  
**Additional Sites:** IUC-TOULOUSE  
**Additional Site Codes:** IUC-TOULOUSE  
**Patient Position:** HFS

**Created in Version:** Accuray Precision 3.3.1.3 [2]  
**Last Modified by:** STADLER, Marine  
**Saved Deliverable by:** MODESTO, Anouchka  
**Plan Key:** 420634/2  
**Delivery Signature:** KQEXp cdeQA

### Auto Generated Results:

Not Applicable

## Report Snapshots



Imaging Series

Primary Series:Yes

Modality:CT

Scan Date:08 Jun 2023, 04:47:57 PM (hr:min:sec)

Scanner Model:Optima CT580

Patient Position:HFS

Number of Slices:203

| Series UID/Description                               | Study UID/Description:                               |
|--|--|
| 1.2.840.113619.2.278.3.3523880722.264.1686200399.497 | 1.2.840.113619.2.278.3.3523880722.264.1686200399.392 |
| PELVIS   | TDM therapie APC                                     |

|                    | X        | Y        | Z       |
|--------------------|----------|----------|---------|
| Origin (IECp, mm)  | -250.000 | -265.000 | 250.000 |
| Size (IECp, voxel) | 512      | 203      | 512     |
| Spacing (IECp, mm) | 0.98     | 2.50     | 0.98    |